# **Smart Fitness Scale**

Model: ESF00+ Series

Questions or Concerns? support.eu@etekcity.com

### Table of Contents

EN

5
8
9
10
11
12
12
13
16
18
18

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### **EN** Specifications

Weight Capacity	180 kg / 28 st / 400 lb
Graduation	0.05 kg / 0.1 lb
Weight Units	kg / st / lb
Platform	ITO-coated glass
Dimensions	375 x 300 x 24 mm / 14.8 x 11.8 x 1.0 in
Battery Type	1000mAh Li-ion battery
Charging Time	3 hours
Input	DC 5V, 500mA
Auto-Off	30 seconds

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# **Safety Information**

#### Please read and follow all instructions and safety guidelines in this manual.

### Medical

- This scale should not be used to diagnose or treat any medical condition.
- Scale data should not be used as a substitute for medical advice from a healthcare professional.
- This scale uses a harmless electrical current to measure body fat. Consult a healthcare professional before using the smart scale if you use medical electronic equipment such as a pacemaker.

### General Use

- Check the scale before using. **Do not** use the scale if it is damaged in any way.
- The glass platform is slippery when wet.
   Always make sure the platform and your feet are dry before stepping on for measurement.
- Place the scale on a hard, flat surface to avoid tipping.
- To weigh safely and accurately, stand with your feet shoulder-width apart while on the scale. **Do not** stand on the edge of the scale or jump on the scale.

- Avoid standing on the scale for a long period of time, or when not weighing yourself.
- **Do not** overload the scale. Maximum weight capacity is 180 kg / 28 st / 400 lb.
- Keep the scale in a cool area with low humidity. Keep the scale away from heat sources, such as ovens or heaters. Avoid temperature fluctuations. Store the scale in a horizontal position, not vertical.
- **Do not** use while charging. Remove the power cord before using.
- Handle with care. Not for commercial use. Household use only.

# SAVE THESE INSTRUCTIONS

## **Function Diagram**

# **Digital Display**

- A. WiFi Icon
- B. Bluetooth® Icon
- C. Battery Life Indicator
- D. Weight Units
- E. Weight Measurement
- F. Fat/Muscle/Water/Bone/BMR Measurement
- G. BMI Measurement
- H. BMI Tooltips

# VeSync App Setup

**Note:** Due to app updates and improvements, the contents of this manual may change without notice. The VeSync app user interface may appear slightly different.

 To download the VeSync app, scan the QR code, or search "VeSync" in the Apple App Store<sup>®</sup> or Google Play<sup>™</sup> Store.

#### Note:

- This scale uses the VeSync app, not the VeSyncFit app.
- For Android users, you may need to select **Allow** to use VeSync.
- Open the VeSync app. If you already have an account, tap Log In. To create a new account, tap Sign Up.
- 3. Make sure your phone has Bluetooth® turned on.
- **4.** Follow the in-app instructions to add your smart scale.

# **Using the Smart Scale**

#### Note:

- To switch between kilograms (kg), stone (st), and pounds (lb), press the unit switching button on the underside of the scale.
- If you are pregnant, scale measurements will not be accurate.

### Weighing without the VeSync App

- 1. Place the scale on a hard, flat surface. [*Figure 1.1*]
- **2.** Step onto the scale. The scale will automatically turn on and start measuring your weight.
- **3.** Continue standing on the scale until the digits on the screen flash 3 times, showing your final measurement. [*Figure 1.2*]

Figure 1.2

## Maintenance

### **Recalibrating the Scale**

If the scale has been moved or flipped upside down, it must be recalibrated to ensure accurate results.

- 1. Place the scale on a hard, flat surface.
- **2.** Step on the scale with 1 foot until digits appear on the display. [*Figure 3.1*] Then, step off.
- **3.** The scale will show " [ " while calibrating, then " []]]" when complete. [*Figure 3.2*]

Figure 3.1

### **Cleaning the Scale**

- Use a slightly damp cloth to clean the surface of the scale, followed by a soft, dry cloth to wipe away any remaining moisture.
- **Do not** use chemical or abrasive cleaning agents to clean the scale.

### Charging the Battery

The display will show "  $\lfloor \frac{1}{2}$  " when the battery needs to be charged.

- 1. Insert the USB charging cable into the USB charging socket. [*Figure 4.1*]
- 2. Plug the USB charging cable into a DC 5V adapter and plug into an outlet. You can also plug the cable into a powered USB outlet. [*Figure 4.1*]

**3.** The bars in the battery symbol ( 🛲 ) will continuously fill as the scale charges. When charged, the symbol will be filled.

#### Note:

- Charging time is about 3 hours.
- For best results, fully charge your scale before using.

Figure 4.1

### **Display Messages**

### Troubleshooting

#### 1. The scale doesn't turn on.

• The scale may need to be charged. Plug it in and charge for 3 hours.

#### 2. Nothing happens when I press the unit button.

• The scale may need to be reset. Use a safety pin or paper clip to press the "Reset" button on the underside of the scale.

# 3. Information flashes on the screen and then disappears.

• The battery is low and may need to be charged.

#### 4. The scale doesn't connect with the VeSync app.

- Make sure your phone's operating system is running on iOS<sup>™</sup> 9.0+ or Android<sup>™</sup> 4.3+.
- Close and relaunch the VeSync app. Make sure you're using the latest version.
- Turn on Bluetooth on your phone (for Android users, make sure your Location is turned on as well).
- Make sure your scale and phone are within 10 m / 30 ft of each other.
- Make sure the scale is not currently connected to any other phone or app.
- Try turning off any VPN apps.
- The scale may need to be reset. Use a safety pin or paper clip to press the "Reset" button on the underside of the scale.

- 5. My scale will only measure my weight, or only my weight and BMI.
  - Set up your smart scale in the VeSync app. Your BMI, Fat, Muscle, Water, Bone, and BMR measurements will show on the display, and more results will show on the app.
  - Make sure Bluetooth is turned on.
  - If you're wearing shoes or socks, the scale cannot measure health metrics such as body fat. Step on the scale with bare feet. [*Figure 5.1*]

**Note:** Resetting the scale will fix many issues.

- A. To clear user data from the scale, press and hold the unit button on the underside of the scale for 10 seconds.
- **B.** To hard reset the scale, use a safety pin or paper clip to press the "Reset" button on the underside of the scale.

If you still need help, please contact **Customer Support** (see page 18).

Conductive Areas

Figure 5.1

## **Warranty Information**

### **Terms & Policy**

Etekcity Corporation warrants all products to be of the highest quality in material, craftsmanship, and service for 2 years, effective from the date of purchase to the end of the warranty period. Warranty lengths may vary between product categories.

Should you encounter any issues or have any questions regarding your new product, feel free to contact our helpful Customer Support Team.

# **Customer Support**

#### **Etekcity Corporation**

1202 N. Miller St., Suite A Anaheim, CA 92806 USA

Email: support.eu@etekcity.com

\*Please have your invoice and order ID ready before contacting Customer Support.

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